

## Self-Seal Autoclave & Sterilisation Pouches

### For Steam Autoclaves & Ethylene Oxide Sterilisers

Excelsior Code: L850/\*E

#### Features

##### **Autoclave Pouches**

Excelsior autoclave pouches have been designed with a special, easy to use, self seal strip, to allow for the quick and efficient autoclaving of individual items.

- pouches are self-seal: peel back, fold over and seal, no separate heat sealer required
- manufactured from 2 ply laminate, transparent film, compliant with EN868-5 and medical grade paper compliant with N868-3
- autoclavable to 121°C and 134°C
- indicators for steam and ethylene sterilisation
- designed with off-set webs for easy opening and tack seals to prevent curling and reduce external Contamination
- the minimum width of any seal is 6mm, tested in accordance to the requirements of Annex C EN868-5
- tamper evident
- alcohol-free
- latex free

Excelsior autoclave pouches are manufactured in a Clean Room facility, which is specified to operate at ISO Class 7 (FS 209E Class 10,000) classification, in accordance with International Standard ISO 14644-1:1999.

#### Range of Sizes Available

PRODUCT CODE	SIZE	POUCHES PER BOX	BOXES PER CASE
L850/AE	90mm x 230mm	200	20
L850/BE	135mm x 255mm	200	15
L850/CE	190mm x 330mm	200	10
L850/DE	300mm x 460mm	200	5



For additional product information:  
Please visit us at  
[www.excelsiorscientific.com](http://www.excelsiorscientific.com)  
Email us at [sales@excelsiorscientific.com](mailto:sales@excelsiorscientific.com)